

REUSABLE STANDARD CAPACITY Oval Cup

A

Old Product Codes	Description	Jaw O.D (mm)	Working Length (cm)	Minimum Working Channel (mm)	Color Code	Packaging
BF-RA-2460	Standard Oval Cup	2.4	600	φ 2.8	●	2/Box
BF-RA-2412	Standard Oval Cup	2.4	1200	φ 2.8	●	2/Box
BF-RA-2416	Standard Oval Cup	2.4	1600	φ 2.8	●	2/Box
BF-RA-2418	Standard Oval Cup	2.4	1800	φ 2.8	●	2/Box
BF-RA-2420	Standard Oval Cup	2.4	2000	φ 2.8	●	2/Box
BF-RA-2423	Standard Oval Cup	2.4	2300	φ 2.8	●	2/Box
BF-RA-2427	Standard Oval Cup	2.4	2700	φ 2.8	●	2/Box
BF-RA-1860	Standard Oval Cup	1.8	600	φ 2.0	●	2/Box
BF-RA-1812	Standard Oval Cup	1.8	1200	φ 2.0	●	2/Box
BF-RA-1816	Standard Oval Cup	1.8	1600	φ 2.0	●	2/Box
BF-RA-1818	Standard Oval Cup	1.8	1800	φ 2.0	●	2/Box

REUSABLE STANDARD CAPACITY Oval Cup w/Needle

B

Old Product Codes	Description	Jaw O.D (mm)	Working Length (cm)	Minimum Working Channel (mm)	Color Code	Packaging
BF-RB-2460	Standard Oval Cup w/Needle	2.4	600	φ 2.8	●	2/Box
BF-RB-2412	Standard Oval Cup w/Needle	2.4	1200	φ 2.8	●	2/Box
BF-RB-2416	Standard Oval Cup w/Needle	2.4	1600	φ 2.8	●	2/Box
BF-RB-2418	Standard Oval Cup w/Needle	2.4	1800	φ 2.8	●	2/Box
BF-RB-2420	Standard Oval Cup w/Needle	2.4	2000	φ 2.8	●	2/Box
BF-RB-2423	Standard Oval Cup w/Needle	2.4	2300	φ 2.8	●	2/Box
BF-RB-2427	Standard Oval Cup w/Needle	2.4	2700	φ 2.8	●	2/Box

REUSABLE STANDARD CAPACITY Alligator Cup

C

Old Product Codes	Description	Jaw O.D (mm)	Working Length (cm)	Minimum Working Channel (mm)	Color Code	Packaging
BF-RC-2460	Standard Alligator Cup	2.4	600	φ 2.8	●	2/Box
BF-RC-2412	Standard Alligator Cup	2.4	1200	φ 2.8	●	2/Box
BF-RC-2416	Standard Alligator Cup	2.4	1600	φ 2.8	●	2/Box
BF-RC-2418	Standard Alligator Cup	2.4	1800	φ 2.8	●	2/Box
BF-RC-2420	Standard Alligator Cup	2.4	2000	φ 2.8	●	2/Box
BF-RC-2423	Standard Alligator Cup	2.4	2300	φ 2.8	●	2/Box
BF-RC-2427	Standard Alligator Cup	2.4	2700	φ 2.8	●	2/Box
BF-RC-1860	Standard Alligator Cup	1.8	600	φ 2.0	●	2/Box
BF-RC-1812	Standard Alligator Cup	1.8	1200	φ 2.0	●	2/Box
BF-RC-1816	Standard Alligator Cup	1.8	1600	φ 2.0	●	2/Box
BF-RC-1818	Standard Alligator Cup	1.8	1800	φ 2.0	●	2/Box

REUSABLE STANDARD CAPACITY Alligator Cup w/Needle

D

Old Product Codes	Description	Jaw O.D (mm)	Working Length (cm)	Minimum Working Channel (mm)	Color Code	Packaging
BF-RD-2460	Standard Alligator Cup w/ Needle	2.4	600	φ 2.8	●	2/Box
BF-RD-2412	Standard Alligator Cup w/ Needle	2.4	1200	φ 2.8	●	2/Box
BF-RD-2416	Standard Alligator Cup w/ Needle	2.4	1600	φ 2.8	●	2/Box
BF-RD-2418	Standard Alligator Cup w/ Needle	2.4	1800	φ 2.8	●	2/Box
BF-RD-2420	Standard Alligator Cup w/ Needle	2.4	2000	φ 2.8	●	2/Box
BF-RD-2423	Standard Alligator Cup w/ Needle	2.4	2300	φ 2.8	●	2/Box
BF-RD-2427	Standard Alligator Cup w/ Needle	2.4	2700	φ 2.8	●	2/Box